
**MANUAL  
HANDLEIDING  
MODE D'EMPLOI  
ANLEITUNG**

## peqMETER 1.14

Input pH/mV: pH, mV

Input EC: Conductivity

Input °C: Temperature (Pt1000)

## Index

General informations.....	1
Installation of the instrument.....	2
pH/mV/temperature measurement.....	2
Conductivity/TDS/Salinity measurement.....	3
Maintenance of electrodes.....	4
Good Laboratory Practice.....	5
Data Storage in the internal memory.....	5
 Algemene informatie.....	 7
Instelling van het toestel.....	8
pH/mV/temperatuur-meting.....	8
Geleidbaarheid/TDS/Saliniteitsmeting.....	9
Onderhoud van de elektroden.....	10
Goede Laboratorium Praktijk.....	11
Data-opslag in het intern geheugen.....	11
 Informations générales.....	 13
Installation de l'appareil.....	14
Mesure de pH/mV/Température.....	14
Mesure de Conductivité/TDS/Salinité.....	15
Entretien des électrodes.....	16
Bonnes Pratiques de Laboratoire.....	17
Collection de données dans la mémoire interne.....	17
 Allgemeine Informationen.....	 19
Installation des Gerätes.....	20
pH/mV/Temperatur-Messung.....	20
Leitfähigkeit/TDS/Salinitätsmessung.....	21
Wartung der Elektroden.....	22
Gute Laborpraxis.....	23
Datenspeicherung im internen Speicher.....	23
 NIST pH buffers.....	 24
Standard pH buffers.....	24
Conductivity standards.....	24
Redox standards.....	24

## General informations

---

### Introduction

This instrument is manufactured with the latest technology and needs no particular maintenance. **PEQLAB** certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 0° and 40°C. Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

### Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. **PEQLAB** will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from **PEQLAB** before returning any product for warranty repair on a freight prepaid basis!

**PEQLAB** is not liable for consequential damages arising out of the use or handling of its products.

### Servicing

In the event of this instrument being returned for servicing, the owner is requested **NOT** to send the following items unless they are suspect:

*Manual*

*Cables*

*Accessories*

If serious malfunctioning occurs, stop using the unit immediately and consult your local **Consort** dealer.

### AC adaptor

Connect the mains adapter to the DC socket. Do not hold the adapter by wet hand.

### Keyboard

**MODE** = Selects all modes or escapes from error traps, calibration procedures, etc. by returning to the original mode.

**CAL** = Starts or proceeds a calibration or a function.

**↑↓** = Button for entering a value or for selecting a function.

**ON/OFF** = Switches the instrument on or off.

### Contact and support

PEQLAB Biotechnologie GmbH

Carl-Thiersch-Str. 2B

91052 Erlangen, Germany

Fon: +49 9131 61070-20

Fax: +49 9131 61070-99

info@peqlab.com

www.peqlab.com

### Inputs

The measuring electrodes should be connected to the corresponding coaxial pH/mV or EC connectors. Automatic temperature compensation and temperature measurements are possible by plugging a Pt1000 temperature probe into the °C terminals. Without Pt1000, the manual temperature compensation is automatically switched on.

- *Never immerse a pH electrode and a dissolved oxygen electrode together in the same solution!*
- *A blinking decimal point warns you for unstable measurements. Wait to read the display! Stirring the solution during the measurements promotes the homogeneity and is obligatory!*
- *The instrument will refuse automatic calibration when the electrode is unstable. Insufficient stirring or a worn electrode may be the cause.*

### Set-up

1. Select [SETTINGS] by pressing **MODE**.
2. Follow the instructions on the screen to adjust language, contrast and automatic power-off timer or to erase the memory. Press **CAL** to confirm or **MODE** to exit.

## pH/mV/temperature measurement

---

### pH measurement

1. Select the desired mode (pH) by pressing **MODE**. The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. The instrument permits to choose between the buffers in memory (1.68, 2.00, 4.00, 4.01, 6.87, 7.00, 9.18, 9.21, 10.01, 12.00, 12.45). Select the proper values and press **CAL**. The unused buffers should be switched off.
3. Rinse the electrodes with distilled water and immerse them in the first buffer solution. Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
4. After rinsing the electrodes with distilled water, immerse them in the samples and read the display.
5. Rinse the electrodes always with distilled water after use and store them in a 3...4 M KCl solution.

### mV measurement

1. Select the desired mode (mV) by pressing **MODE**. The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
  2. Immerse the electrodes in a standard solution of known potential. Select [CALIBRATE] and press **CAL**. Adjust to the proper value and press **CAL**.
- *Select [RESET] and press CAL to reset the calibration.*

### Temperature measurement

1. Select the desired mode (°C) by pressing **MODE**. Should you want to recalibrate, press **CAL**.
  2. Immerse the Pt1000 in a solution of known temperature. Select [CALIBRATE] and press **CAL**. Adjust to the proper value and press **CAL**.
- *Select [RESET] and press CAL to reset the calibration.*
  - *Without Pt1000, press CAL, adjust the manual temperature compensation and proceed by pressing CAL.*

## **Conductivity/TDS/Salinity measurement**

---

### **Conductivity measurement**

1. Select the desired mode (S/cm) by pressing **MODE**. The display will immediately show the measured value according to the previous calibration. Should you want to recalibrate, press **CAL**.
2. The instrument permits to choose between the standards in memory (0.01, 0.1, 1 M KCl). Select the proper values and press **CAL**. The unused standards should be switched off.
3. After rinsing the electrode several times with the first standard solution, immerse it in that standard. The solution temperature is not so critical but should lie between 0°C and 40°C. When no Pt1000 is used, do not forget to compensate manually first! Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
4. Rinse the electrode several times with the sample, immerse it in that solution and read the display.
5. Rinse the electrode always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

### **TDS measurement**

1. Select the desired mode (mg/l)(TDS) by pressing **MODE**. The display will immediately show the measured value.

### **Salinity measurement**

1. Select the desired mode (SAL) by pressing **MODE**. The display will immediately show the measured value.

### pH electrodes

A pH electrode is active and stable only after wetting! For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution before usage. When not in use for shorter periods the electrode should always be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the electrode is not used for a longer period (e.g. one month), please check the closing cap regularly and refill with 3 - 4 M KCL solution if necessary. Always plug the closing cap on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

**Avoid a low pressure inside the electrode!** Therefore always remove the closure from the refilling aperture during the measurements as well as during the calibration. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a cleaning solution overnight and then be cleaned with distilled water before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

### Metal electrodes

Metal electrodes (Pt, Ag, Au) are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

- *Silver electrodes are immersed in a concentrated ammonia solution during one hour.*
- *Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.*

### Conductivity electrodes

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

## ***Good Laboratory Practice***

---

1. Select the desired range by pressing **MODE**.
2. Press **CAL**, select GLP and press **CAL** to display a complete calibration report.

## ***Data Storage in the internal memory***

---

### **Storing data**

1. Press **↓** to store the displayed value in memory. The display shows a next storage address.
1. Select the desired address and press **CAL** to enter the measured value in memory.

### **Recalling data**

1. Press **↑** to recall the stored data. The display shows a list of stored values.
2. Select the desired address to read the stored values.
3. Press **MODE** to return to the normal measurements.

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